

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SCN10A Polyclonal Antibody

多克隆抗体

产品货号: PA10220

第三版

描述:The protein encoded by this gene is a tetrodotoxin-resistant voltage-gated sodium channel alpha subunit. The properties of the channel formed by the encoded transmembrane protein can be altered by interaction with different beta subunits. This protein may be involved in the onset of pain associated with peripheral neuropathy.

抗原:Synthetic peptide of human SCN10A

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in the dorsal root ganglia and sciatic nerve.