

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ITPR3 Polyclonal Antibody

多克隆抗体

产品货号: PA10225

第三版

**描述:** This gene encodes a receptor for inositol 1,4,5-trisphosphate, a second messenger that mediates the release of intracellular calcium. The receptor contains a calcium channel at the C-terminus and the ligand-binding site at the N-terminus. Knockout studies in mice suggest that type 2 and type 3 inositol 1,4,5-trisphosphate receptors play a key role in exocrine secretion underlying energy metabolism and growth.

**抗原:** Synthetic peptide of human ITPR3

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Expressed in intestinal crypt and villus epithelial cells.