

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human KIF22 Polyclonal Antibody

多克隆抗体

产品货号: PA10276

第三版

描述:The protein encoded by this gene is a member of the kinesin-like protein family. The family members are microtubule-dependent molecular motors that transport organelles within cells and move chromosomes during cell division. The C-terminal half of this protein has been shown to bind DNA. Studies with the Xenopus homolog suggests its essential role in metaphase chromosome alignment and maintenance. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

抗原:Synthetic peptide of human KIF22

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in bone, cartilage, joint capsule, ligament, skin, and primary cultured chondrocytes.1 Publication