

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PTPMT1 Polyclonal Antibody

多克隆抗体

产品货号: PA10300

第三版

描述:PTPMT1 (protein tyrosine phosphatase mitochondrial 1), also known as MOSP or PLIP (phosphoinositide lipid phosphatase) and previously known as DUSP23, is a widely expressed PTP membrane protein with high expression levels in pancreatic cells. PTPMT1 exclusively localizes to the matrix face of the inner membrane of the mitochondrion. It is responsible for dephosphorylating mitochondrial proteins and therefore plays a significant role in the production of ATP and secretion of insulin.

抗原:Recombinant protein of human PTPMT1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: