

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CERS3 Polyclonal Antibody

多克隆抗体

产品货号: PA10304

第三版

**描述:** LASS3 is a 383 amino acid membrane protein almost exclusively expressed in testis, suggesting that LASS3 plays an important role in proper testis function. It is also weakly expressed in skin. A transcriptional variant of LASS3 cDNA exists and can result in the production of a 419 amino acid protein (LASS3-long). LASS3 overproduction raises the level of several ceramide species, including C18:0- and C24:0-ceramides.

**抗原:** Synthetic peptide of human CERS3

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Expressed in the epidermis, where it localizes at the interface between the stratum granulosum and the stratum corneum (at protein level).<sup>1</sup> Publication