

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat PTPN5 Polyclonal Antibody

多克隆抗体

产品货号: PA10311

---

第三版

**描述:**Protein tyrosine phosphatase non-receptor type 5 is an enzyme that in humans is encoded by the gene. May regulate the activity of several effector molecules involved in synaptic plasticity and neuronal cell survival, including MAPKs, Src family kinases and NMDA receptors.

**抗原:**Recombinant protein of human PTPN5

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**