

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MAP1LC3C Polyclonal Antibody

多克隆抗体

产品货号: PA10315

第三版

描述:Autophagy is a highly regulated bulk degradation process that plays an important role in cellular maintenance and development. MAP1LC3C is an ortholog of the yeast autophagosome protein Atg8. Ubiquitin-like modifier that plays a crucial role in antibacterial autophagy (xenophagy) through the selective binding of CALCOCO2. Recruits all ATG8 family members to infecting bacteria such as *S. Typhimurium*.

抗原:Synthetic peptide of human MAP1LC3C

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Most abundant in placenta, lung and ovary.