

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MAGEC1 Polyclonal Antibody

多克隆抗体

产品货号: PA10368

第三版

描述: This gene is a member of the melanoma antigen gene (MAGE) family. The proteins of this family are tumor-specific antigens that can be recognized by autologous cytolytic T lymphocytes. This protein contains a large number of unique short repetitive sequences in front of the MAGE-homologous sequence, and therefore is about 800 aa longer than the other MAGE proteins.

抗原: Synthetic peptide of human MAGEC1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: