

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat NDRG3 Polyclonal Antibody

多克隆抗体

产品货号: PA10397

---

第三版

**描述:**The N-myc downstream regulated gene (NDRG) family is comprised of four members, NDRG1 (also designated Drg1, RTP, rit42, Cap43 and Ndr1), NDRG2, NDRG3 and NDRG4, which share 57-65% homology. The gene encoding NDRG3 maps to human chromosome 20q11 and is predominantly expressed in testis, prostate and ovary, which suggests it may play a role in spermatogenesis.

**抗原:**Synthetic peptide of human NDRG3

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**Ubiquitous. Highly expressed in brain.