

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ERRF1 Polyclonal Antibody

多克隆抗体

产品货号: PA10426

第三版

描述:ERRF1 is a cytoplasmic protein whose expression is upregulated with cell growth (Wick et al., 1995 [PubMed 7641805]). It shares significant homology with the protein product of rat gene-33, which is induced during cell stress and mediates cell signaling

抗原:Synthetic peptide of human ERRF1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: