

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ROBO4 Polyclonal Antibody

多克隆抗体

产品货号: PA10440

第三版

描述:Robo4 (roundabout homolog 4 or magic roundabout) is a receptor for Slit proteins which implicate in angiogenesis and vascular patterning especially through facilitating inhibition of primary endothelial cell migration. It is a fairly recently identified member of Robo family and is shorter than other Robo members at both intracellular as well as extracellular domains. Robo4 plays important roles in several aspects of vascular development such as the guidance of endothelial cell migration, cell cycle regulation, inhibition of pathologic angiogenesis and endothelial hyperpermeability etc. Employing mouse and zebrafish models, Slit2/Robo4 interactions have been shown to maintain integrity of vascular network and its barrier function by inhibiting cytokine-mediated vasculogenesis with enhanced permeability, and that Slit2-Robo4-paxillin-GIT1 network inhibits neovascularization and vascular leakage. Moreover, single nucleotide polymorphism of Robo4 has been documented to be associated with autism suggesting its role in development and function of human brain.

抗原:Recombinant protein of human ROBO4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Specifically expressed in endothelial cells. Expressed at sites of angiogenesis in different tumor types.1 Publication