

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ADRB3 Polyclonal Antibody

多克隆抗体

产品货号: PA10530

第三版

描述:The protein encoded by this gene belongs to the family of beta adrenergic receptors, which Human, Mouse induce catecholamine-induced activation of adenylate cyclase through the action of G proteins. This receptor is located mainly in the adipose tissue and is involved in the regulation of lipolysis and thermogenesis.

抗原:Synthetic peptide of human ADRB3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed mainly in adipose tissues.