

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SLC25A26 Polyclonal Antibody

多克隆抗体

产品货号: PA10534

第三版

描述: Mitochondrial carriers, including SLC25A26, are a family of transport proteins found mostly in the inner membranes of mitochondria. They shuttle metabolites and cofactors through the mitochondrial membrane

抗原: Recombinant protein of human SLC25A26

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed. Highly expressed in testis, with moderate expression in brain, heart, kidney, lung, skeletal muscle, pancreas, small intestine and liver, and low expression in spleen.