
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CLEC16A Polyclonal Antibody

多克隆抗体

产品货号: PA10548

第三版

描述: This gene encodes a member of the C-type lectin domain containing family. Single nucleotide polymorphisms in introns of this gene have been associated with diabetes mellitus, multiple sclerosis and rheumatoid arthritis. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原: Synthetic peptide of human CLEC16A

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Almost exclusively expressed in immune cells, including dendritic cells, B-lymphocytes and natural killer cells.1 Publication