

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human C12orf75 Polyclonal Antibody

多克隆抗体

产品货号: PA10582

第三版

**描述:**C12orf75 (chromosome 12 open reading frame 75), also known as AGD3 or OCC1, is a protein that in humans is encoded by the C12orf75 gene. C12orf75 is high expression in placenta, skeletal muscle, kidney and pancreas tissues.

**抗原:**Synthetic peptide of human C12orf75

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**