

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ST18 Polyclonal Antibody

多克隆抗体

产品货号: PA10600

第三版

**描述:**ST18 (suppression of tumorigenicity 18), also known as ZNF387 ((zinc finger protein 387), is a 1,047 amino acid nuclear protein that is expressed at low levels in heart, liver, kidney, skeletal muscle, pancreas, testis, ovary and prostate. Containing six C2HC-type zinc fingers and an SMC domain, ST18 is a repressor protein that binds to DNA sequences containing a specific bipartite element. ST18 is considered a transcription factor that represses basal transcription activity from target promoters and Inhibits colony formation in cultured breast cancer cells.

**抗原:**Recombinant protein of human ST18

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Detected at low levels in heart, liver, kidney, skeletal muscle, pancreas, testis, ovary and prostate. Detected at even lower levels in mammary epithelial cells and breast cancer cells.