

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TP53AIP1 Polyclonal Antibody

多克隆抗体

产品货号: PA10625

第三版

描述: This gene is specifically expressed in the thymus, and encodes a protein that is localized to the mitochondrion. The expression of this gene is inducible by p53, and it is thought to play an important role in mediating p53-dependent apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

抗原: Synthetic peptide of human TP53AIP1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Only found to be expressed in thymus.