

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SULT2B1 Polyclonal Antibody

多克隆抗体

产品货号: PA10637

第三版

描述: Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. This gene sulfates dehydroepiandrosterone but not 4-nitrophenol, a typical substrate for the phenol and estrogen sulfotransferase subfamilies. Two alternatively spliced variants that encode different isoforms have been described.

抗原: Recombinant protein of human SULT2B1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed highly in placenta, prostate and trachea and lower expression in the small intestine and lung.