

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat STX7 Polyclonal Antibody

多克隆抗体

产品货号: PA10649

第三版

描述:The Syntaxin 7 antibody from Proteintech is a rabbit polyclonal antibody to a recombinant protein of human Syntaxin 7. This antibody recognizes human,mouse, rat antigen. The Syntaxin 7 antibody has been validated for the following applications: ELISA, WB, IHC, IP analysis.

抗原:Recombinant protein of human STX7

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highest expression is found in placenta followed by heart, skeletal muscle, kidney and brain. Low expression is found in pancreas, lung and liver.