

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human STX10 Polyclonal Antibody

多克隆抗体

产品货号: PA10660

第三版

**描述:**This gene belongs to the syntaxin family and encodes a soluble N-ethylmaleimide sensitive factor attachment protein receptor (SNARE). The encoded protein is involved in docking and fusion events at the Golgi apparatus. Alternative splicing results in multiple transcript variants.

**抗原:**Recombinant protein of human STX10

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed at high levels in heart, skeletal muscle and pancreas.