

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR92 Polyclonal Antibody

多克隆抗体

产品货号: PA10676

第三版

描述:This gene encodes a member of the rhodopsin class of G protein-coupled transmembrane receptors. This protein transmits extracellular signals from lysophosphatidic acid to cells through heterotrimeric G proteins and mediates numerous cellular processes. Many G protein receptors serve as targets for pharmaceutical drugs. Transcript variants of this gene have been described.

抗原:Synthetic peptide of human LPAR5

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Not expressed in frontal cortex, basal forebrain, caudate putamen, thalamus, or hippocampus.