

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TH Polyclonal Antibody

多克隆抗体

产品货号: PA10689

第三版

描述:The protein encoded by this gene is involved in the conversion of tyrosine to dopamine. It is the rate-limiting enzyme in the synthesis of catecholamines, hence plays a key role in the physiology of adrenergic neurons. Mutations in this gene have been associated with autosomal recessive Segawa syndrome. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

抗原:Synthetic peptide of human TH

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Mainly expressed in the brain and adrenal glands.