

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse RAB39 Polyclonal Antibody

多克隆抗体

产品货号: PA10733

第三版

描述: This gene encodes a small GTP-binding protein that belongs to the largest family within the Ras superfamily. Rab39 is a Golgi-associated protein involved in endocytosis of HeLa cells. It contains conserved motifs referred to phosphate/guanosine binding and a microbody C-terminal targeting signal. Rab39 is mostly detected in epithelial cell lines and expressed ubiquitously in human tissues. It can facilitate endocytosis in HeLa cells when expressed transiently and stably.

抗原: Synthetic peptide of human RAB39

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: