

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RAB41 Polyclonal Antibody

多克隆抗体

产品货号: PA10736

第三版

描述: This gene encodes a small GTP-binding protein that belongs to the largest family within the Ras superfamily. These proteins function as regulators of membrane trafficking. They cycle between inactive GDP-bound and activated GTP-bound states, which is controlled by GTP hydrolysis-activating proteins (GAPs). This family member can be activated by the GAP protein RN-Tre, and it is localized to the Golgi complex.

抗原: Synthetic peptide of human RAB41

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed in brain, testis, lung, heart, ovary, colon, kidney, uterus and spleen but not in liver.1 Publication