

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SPAG1 Polyclonal Antibody

多克隆抗体

产品货号: PA10745

第三版

描述:The correlation of anti-sperm antibodies with cases of unexplained infertility implicates a role for these antibodies in blocking fertilization. Improved diagnosis and treatment of immunologic infertility, as well as identification of proteins for targeted contraception, are dependent on the identification and characterization of relevant sperm antigens. The protein expressed by this gene is recognized by anti-sperm agglutinating antibodies from an infertile woman. Furthermore, immunization of female rats with the recombinant human protein reduced fertility. This protein localizes to the plasma membrane of germ cells in the testis and to the post-acrosomal plasma membrane of mature spermatozoa. Recombinant polypeptide binds GTP and exhibits GTPase activity. Thus, this protein may regulate GTP signal transduction pathways involved in spermatogenesis and fertilization. Two transcript variants of this gene encode the same protein.

抗原:Synthetic peptide of human SPAG1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Present in most tissues, including lung, with the strongest expression in brain, colon, kidney, and testis. In sperm and testis, detected in particular in pachytene primary spermatocytes. Up-regulated in pancreatic tumor tissues and not in normal pancreatic tissue.