
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RAB5A Polyclonal Antibody

多克隆抗体

产品货号: PA10761

第三版

描述: Ras-related protein Rab-5A is a protein that in humans is encoded by the RAB5A gene. It is rate-limiting for homotypic early endosome fusion. It required for the fusion of plasma membranes and early endosomes, as well as contributes to the regulation of filopodia extension.

抗原: Synthetic peptide of human RAB5A

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: