

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RAB3D Polyclonal Antibody

多克隆抗体

产品货号: PA10767

第三版

描述:The RAB3D cDNA encodes a 24-kD protein as determined by in vitro translation. The deduced 219-amino acid sequence shows highest homology to the RAB3 subfamily of low molecular mass GTP-binding proteins. Northern blot analysis of various human hematopoietic cell lines showed a single transcript of 4.1 kb present at variable levels in each cell line, with highest expression in granulocytic lines. RAB3D was also highly expressed in YJ and HL-60 cells in response to retinoic acid treatment.

抗原:Synthetic peptide of human Rab3D

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in granulocytes of peripheral blood. Constitutively expressed at low levels in all hematopoietic cell lines investigated.