

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PTPN20 Polyclonal Antibody

多克隆抗体

产品货号: PA10778

第三版

描述:The product of this gene belongs to the family of classical tyrosine-specific protein tyrosine phosphatases. Many protein tyrosine phosphatases have been shown to regulate fundamental cellular processes. The encoded protein appears to be targeted to sites of actin polymerization. A pseudogene of this gene has been defined on chromosome 10. Alternative splicing results in multiple transcript variants.

抗原:Synthetic peptide of human PTPN20A/PTPN20B

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Present in many cell lines (at protein level). Widely expressed.