

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RNF126 Polyclonal Antibody

多克隆抗体

产品货号: PA10793

第三版

描述:The protein encoded by this gene contains a RING finger domain, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions. RNF126 interacts with p21 and RNF126 over-expression increased p21 protein ubiquitination in an E3 ligase activity-dependent manner. RNF126 knockdown induced loss of cell viability in MDA-MB-231 and PC-3 can be partially rescued by depletion of p21. RNF126 stable knockdown in PC3 inhibited tumor growth in SCID mice.

抗原:Synthetic peptide of human RNF126

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in liver and testis.