

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse PPP1R13B Polyclonal Antibody

多克隆抗体

产品货号: PA10839

第三版

描述:The protein plays a central role in regulation of apoptosis via its interaction with p53/TP53.Regulates TP53 by enhancing the DNA binding and transactivation function of TP53 on the promoters of proapoptotic genes in vivo.It contains four ankyrin repeats and an SH3 domain involved in protein-protein interactions. ASPP proteins are required for the induction of apoptosis by p53-family proteins.

抗原:Synthetic peptide of human ASPP1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Reduced expression in breast carcinomas expressing a wild-type TP53 protein.