

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRIM34 Polyclonal Antibody

多克隆抗体

产品货号: PA10870

第三版

描述:The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, B-box type 1 and B-box type 2 domain, and a coiled-coil region. Expression of this gene is up-regulated by interferon. This gene is mapped to chromosome 11p15, where it resides within a TRIM gene cluster. Alternative splicing results in multiple transcript variants.

抗原:Synthetic peptide of human TRIM34

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Isoform 1 is the most abundant form. It is highly expressed in the placenta, spleen, colon and peripheral blood leukocytes.