

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RRAS Polyclonal Antibody

多克隆抗体

产品货号: PA10883

第三版

描述:Ras-related protein R-Ras is a protein that in humans is encoded by the RRAS gene. RRAS has been shown to interact with RASSF5, NCK1, Bcl-2, ARAF and RALGDS. Regulates the organization of the actin cytoskeleton. Interacts with PLCE1. Interacts (active GTP-bound form preferentially) with RGS14.

抗原:Synthetic peptide of human RRAS

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: