

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse AKAP8L Polyclonal Antibody

多克隆抗体

产品货号: PA10886

第三版

描述:A-kinase anchor protein 8-like is an enzyme that in humans is encoded by the AKAP8L gene. AKAP8L has been shown to interact with Thymopoietin[4][5] and DHX9. Could play a role in constitutive transport element (CTE)-mediated gene expression. Does not seem to be implicated in the binding of regulatory subunit II of PKA. May be involved in nuclear envelope breakdown and chromatin condensation. May regulate the initiation phase of DNA replication when associated with TMPO-beta.

抗原:Synthetic peptide of human AKAP8L

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed. Expressed in the brain cortex (at protein level).