

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SAMD4A Polyclonal Antibody

多克隆抗体

产品货号: PA10899

第三版

描述: Sterile alpha motifs (SAMs) in proteins such as SAMD4A are part of an RNA-binding domain that functions as a posttranscriptional regulator by binding to an RNA sequence motif known as the Smaug recognition element, which was named after the *Drosophila* Smaug protein.

抗原: Synthetic peptide of human SAMD4A

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: