

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SAMD9 Polyclonal Antibody

多克隆抗体

产品货号: PA10902

第三版

描述: This gene encodes a sterile alpha motif domain-containing protein. The encoded protein localizes to the cytoplasm and may play a role in regulating cell proliferation and apoptosis. Mutations in this gene are the cause of normophosphatemic familial tumoral calcinosis. Alternate splicing results in multiple transcript variants that encode the same protein.

抗原: Synthetic peptide of human SAMD9

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed. Very low levels in skeletal muscle. Not detected in fetal brain. Down-regulated in aggressive fibromatosis, as well as in breast and colon cancers.