

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD66e Polyclonal Antibody

多克隆抗体

产品货号: PA10935

第三版

**描述:**Members of the CEACAM subfamily, including CEACAM5, belong to the CEA gene family. For general information on the CEA gene family, see CEACAM1. Cell surface glycoprotein that plays a role in cell adhesion and in intracellular signaling. Receptor for E.coli Dr adhesins. Found in adenocarcinomas of endodermally derived digestive system epithelium and fetal colon.

**抗原:**Synthetic peptide of human CEACAM5

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Found in adenocarcinomas of endodermally derived digestive system epithelium and fetal colon.