

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CALCA Polyclonal Antibody

多克隆抗体

产品货号: PA10993

第三版

**描述:** This gene encodes the peptide hormones calcitonin, calcitonin gene-related peptide and katecalcitonin by tissue-specific alternative RNA splicing of the gene transcripts and cleavage of inactive precursor proteins. Calcitonin is involved in calcium regulation and acts to regulate phosphorus metabolism. Calcitonin gene-related peptide functions as a vasodilator while katecalcitonin is a calcium-lowering peptide. Multiple transcript variants encoding different isoforms have been found for this gene.

**抗原:** Synthetic peptide of human CALCA

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**