

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SLC14A1 Polyclonal Antibody

多克隆抗体

产品货号: PA11001

第三版

**描述:**The protein encoded by this gene is a membrane transporter that mediates urea transport in erythrocytes. This gene forms the basis for the Kidd blood group system. Mediates urea transport in erythrocytes By similarity. Low-affinity facilitative urea transporter that allow rapid equilibration between the urinary space and the hyperosmotic interstitium. the rate of urea conduction is increased by hypoosmotic stress.

**抗原:**Synthetic peptide of human SLC14A1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**Detected in erythrocytes (at protein level). Erythrocytes.1 Publication