

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LRP3 Polyclonal Antibody

多克隆抗体

产品货号: PA11003

第三版

描述:Low density lipoprotein receptor-related protein 3 (LRP-3) is a protein that in humans is encoded by theLRP3 gene. Probable receptor, which may be involved in the internalization of lipophilic molecules and/or signal transduction. Its precise role is however unclear, since it does not bind to very low density lipoprotein (VLDL) or to LRPAP1 in vitro.

抗原:Synthetic peptide of human LRP3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed. Highly expressed in skeletal muscle and ovary. Expressed at intermediate level in heart, brain, liver, pancreas, prostate and small intestine. Weakly expressed in testis, colon and leukocyte.