

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BPIFB3 Polyclonal Antibody

多克隆抗体

产品货号: PA11012

第三版

**描述:**BPI fold-containing family B member 3 may have the capacity to recognize and bind specific classes of odorants. May act as a carrier molecule, transporting odorants across the mucus layer to access receptor sites. May serve as a primary defense mechanism by recognizing and removing potentially harmful odorants or pathogenic microorganisms from the mucosa or clearing excess odorant from mucus to enable new odorant stimuli to be received.

**抗原:**Synthetic peptide of human BPIFB3

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$  g/ml.

**特异性:**Detected in nasal septal epithelium.1 Publication