

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC44A2 Polyclonal Antibody

多克隆抗体

产品货号: PA11055

第三版

**描述:** CTL2, also called SLC44A2 (solute carrier family 44 member 2), is a multi-pass membrane protein expressed in cells of the inner ear. CTL2 is a possible candidate for autoimmune hearing loss in humans. Present in supporting cells of the inner ear (at protein level). Only isoform 3 is expressed in inner ear vestibular tissue.

**抗原:** Synthetic peptide of human SLC44A2

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Present in supporting cells of the inner ear (at protein level).