

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SLC5A11 Polyclonal Antibody

多克隆抗体

产品货号: PA11059

第三版

**描述:** Cotransporters, such as SLC5A11, represent a major class of proteins that make use of ion gradients to drive active transport for the cellular accumulation of nutrients, neurotransmitters, osmolytes, and ions Roll et al. (2002) [PubMed 12039040].

**抗原:** Synthetic peptide of human SLC5A11

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Highest expression in heart, skeletal muscle, kidney, liver and placenta. Weaker expression in brain, colon, spleen, lung and peripheral blood leukocytes.