

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD205 Polyclonal Antibody

多克隆抗体

产品货号: PA11061

第三版

**描述:**Lymphocyte antigen 75 is a protein that in humans is encoded by the LY75 gene.To direct captured antigens from the extracellular space to a specialized antigen-processing compartment as an endocytic receptor. Causes reduced proliferation of B-lymphocytes.Selection splicing results in multiple transcript variants encoding fusion products that share sequence identity with each individual gene product.

**抗原:**Synthetic peptide of human LY75

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$  g/ml.

**特异性:**Expressed in spleen, thymus, colon and peripheral blood lymphocytes. Detected in myeloid and B lymphoid cell lines. Isoform 2 and isoform 3 are expressed in malignant Hodgkin lymphoma cells called Hodgkin and Reed-Sternberg (HRS) cells.