

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC5A9 Polyclonal Antibody

多克隆抗体

产品货号: PA11064

第三版

**描述:**SGLT-4 (sodium/glucose cotransporter 4), also known as SLC5A9 (solute carrier family 5 member 9), is a 681 amino acid multi-pass membrane protein that participates in the sodium-dependent transport of D-fructose, D-mannose and D-glucose. A member of the sodium:solute symporter (SSF) family, SGLT-4 is found in small intestine and kidney.

**抗原:**Synthetic peptide of human SLC5A9

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in the small intestine, kidney and liver.1 Publication