

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IRS4 Polyclonal Antibody

多克隆抗体

产品货号: PA11085

第三版

描述: IRS4 encodes the insulin receptor substrate 4, a cytoplasmic protein that contains many potential tyrosine and serine/threonine phosphorylation sites. Tyrosine-phosphorylated IRS4 protein has been shown to associate with cytoplasmic signalling molecules that contain SH2 domains. The IRS4 protein is phosphorylated by the insulin receptor tyrosine kinase upon receptor stimulation

抗原: Synthetic peptide of human IRS4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in myoblasts. Expressed in liver and hepatocellular carcinoma.