

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EDNRB Polyclonal Antibody

多克隆抗体

产品货号: PA11108

第三版

**描述:**The protein encoded by this gene is a G protein-coupled receptor which activates a phosphatidylinositol-calcium second messenger system. Its ligand, endothelin, consists of a family of three potent vasoactive peptides: ET1, ET2, and ET3. Studies suggest that the multigenic disorder, Hirschsprung disease type 2, is due to mutations in the endothelin receptor type B gene. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**抗原:**Synthetic peptide of human EDNRB

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in placental stem villi vessels, but not in cultured placental villi smooth muscle cells.