

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PTGER2 Polyclonal Antibody

多克隆抗体

产品货号: PA11130

第三版

**描述:** This gene encodes a receptor for prostaglandin E2, a metabolite of arachidonic acid which has different biologic activities in a wide range of tissues. Mutations in this gene are associated with aspirin-induced susceptibility to asthma. Prostaglandin E2 receptor (53kDa), also known as EP2, is a prostaglandin receptor for prostaglandin E2. It is encoded by the human gene PTGER2. Receptor for prostaglandin E2 (PGE2). The activity of this receptor is mediated by G(s) proteins that stimulate adenylate cyclase. The subsequent raise in intracellular cAMP is responsible for the relaxing effect of this receptor on smooth muscle.

**抗原:** Synthetic peptide of human PTGER2

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Placenta and lung.