

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HINT2 Polyclonal Antibody

多克隆抗体

产品货号: PA11146

第三版

描述: Histidine triad proteins, such as HINT2, are nucleotide hydrolases and transferases that act on the alpha-phosphate of ribonucleotides. Hydrolase probably involved in steroid biosynthesis. May play a role in apoptosis. Has adenosine phosphoramidase activity. High expression in liver and pancreas. Expression is significantly down-regulated in hepatocellular carcinoma (HCC) patients.

抗原: Synthetic peptide of human HINT2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: High expression in liver and pancreas. Expression is significantly down-regulated in hepatocellular carcinoma (HCC) patients.