

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FPR2 Polyclonal Antibody

多克隆抗体

产品货号: PA11189

第三版

描述:Low affinity receptor for N-formyl-methionyl peptides, which are powerful neutrophils chemotactic factors. Binding of FMLP to the receptor causes activation of neutrophils. This response is mediated via a G-protein that activates a phosphatidylinositol-calcium second messenger system. The activation of LXA4R could result in an anti-inflammatory outcome counteracting the actions of proinflammatory signals such as LTB4

抗原:Synthetic peptide of human FPR2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed abundantly in the lung and neutrophils. Also found in the spleen and testis.