

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SYNPO Polyclonal Antibody

多克隆抗体

产品货号: PA11208

第三版

描述:Synaptopodin is an actin-associated protein that may play a role in actin-based cell shape and motility. The name synaptopodin derives from the protein's associations with postsynaptic densities and dendritic spines and with renal podocytes.

抗原:Synthetic peptide of human SYNPO

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in cerebral cortex.